



Doc No.: RP-RA-03 RE01

## Instructions For use

This set of instructions for use was written with respect to the following non-invasive medical devices:

- Virox 70% Isopropyl Alcohol Wipes 1s (50 x 1s) (Product code: VR3926)
- Virox 70% Isopropyl Alcohol Wipes 1s (100 x 1s) (Product code: VR3927 and VR3927NKF)

The Basic UDI-DI of these non-invasive medical devices is 88863897ViroxV5.

Our intended use of these products is to kill 99.99% of germs. It is designed for disinfection and cleaning of hard surfaces. Not for end-point sterilization. Our indications of use are for disinfection of non-invasive medical devices. Not for end-point sterilisation. These products are designed for disinfection of all surfaces, such as metal, rubber and fabrics, and are meant to be disposed after single use. Avoid acrylic glass, acrylic surfaces and other plastic surfaces when using alcohol-based wipes.

This product is for external use only. It should be stored in a cool place, away from naked flame and direct sunlight. Do not flush. The wipes can be used on general non-invasive medical devices as depicted in Figure 1, such as hospital beds, wheelchairs and thermometers. Make sure that all visible soiling on target surface is removed before disinfection. Wipe surface thoroughly with the impregnated wipe, ensure complete wetting and allow the surface to dry completely. Duration of using wet wipes is dependent on the drying time. Ensure that the antimicrobial efficiency of the product matches its intended use. Not for end-point disinfection of semicritical and critical medical devices.




Figure 1: Examples of non-invasive medical devices that are suitable for wet wipes

As a precaution, users should keep this product out of reach of children. This product is flammable, hence keep away from heat, sparks and open flame. Avoid contact with eyes. Store at room temperature. Wearing of gloves (eg. Nitrile rubber) is recommended. First aid: if in contact with eyes; hold eyes open and rinse slowly and gently with water for 15-20 minutes. If contact lenses are present, remove them immediately and rinse eyes. Call a poison control center for treatment advice.

